

# The latest generation of beds for care at home.

domiflex® 3

.bock®

# New standard for care at home

The new bock domiflex® 3 bed series takes an important step into the future to improve care at home. With our more than 100 years of experience we have further developed our classic care beds domiflex® and domiflex® 2 and combined all advantages in one bed range: domiflex® 3.

Even the standard model has an extended height adjustment range of 35–80 cm, and offers an increased patient weight of up to 155 kg and thus a higher safe working load of 190 kg. At the same time, one can choose between two different colours: the classic Steamed Beech or the modern Sonoma Oak Light.



Can be lowered to 35 cm: domiflex® 3 sets new standards in height adjustment and enables ergonomic and safe entry and exit.



Cable-connected hand control with locking function and numerous adjustment possibilities. The soft buttons provide comfortable haptics and safe handling.



Choice of two colours:  
Here in Steamed Beech.



Coming soon: heavy duty version up to 185 kg patient weight.

# domiflex<sup>®</sup> 3: all advantages at a glance



Extended height adjustment range of 35–80 cm for ergonomic entry and exit and optimal care height



Space saving transport and storage system in combination with light components offer an easy handling



Tool-free connection technology enables a quick and easy assembly and disassembly



Choice of two standard colours: Sonoma Oak Light and Steamed Beech



Comprehensive range of accessories for individual care



Quality made in Germany

# domiflex<sup>®</sup> 3: technical data



## domiflex<sup>®</sup> 3

### Dimensions

External dimensions	103 × 210.7 cm
Lying surface	90 × 200 cm
Special dimensions	Bed extension to 220 cm possible

Height adjustment	35–80 cm
Lifter space	> 15 cm


### Weights

Max. patient weight	155 kg
Safe working load	190 kg
Bed weight	78 kg
Max. weight per component	< 18 kg


### Motoren

Splash-proof according to IPX4, high-quality spiral cable with strain relief and kink protection, SMPS plug transformer with extra-low safety voltage offering safety from the socket.


## Colours



Sonoma Oak Light



Steamed Beech



# More safety, more comfort: domiflex® 3 accessories

## Bed extension to 220 cm



### Set consisting of:

- Frame extension with tool-free connection technology
- 220 cm side rails
- Lying surface extension

## Side rails SR2



### Safe and proverbially transparent

- Can be retrofitted at any time
- Safe and convenient entry and exit
- Modern design



## Mattress compensation



For a comfortable and ergonomic sitting position and reduction of pressure points

## Urine bottle holder



## Drainage bag holder



# More safety, more comfort: domiflex® 3 accessories

## Bed lamp F



- Lighting: LED, 35 V/5 W
- Switchable via hand control

## Hand control with additional functions



This hand control offers additionally Trendelenburg and comfort sitting position

## Angled infusion holder



## Height-adjustable infusion holder



## Infusion holder for lifting pole



# More safety, more comfort: domiflex® 3 accessories

**Grab rail**



**Side rail attachment**




**Side rail bumpers**



**Side rail bumpers**



**Wooden end panel covering**



# domiflex<sup>®</sup> 3: the bock marketing package

## Marketing and sales support

- ⊕ Operation manual as a short version along with each bed
- ⊕ Comprehensive photo and video material
- ⊕ 2-page flyer
- ⊕ bock ACADEMY - product training sessions
- ⊕ Support and content for your marketing activities



**bock videos - YouTube**  
Here you can always find the latest videos about domiflex<sup>®</sup> 3.

You can find the latest news about our company and our products online.  
On our website [www.bock.net](http://www.bock.net) as well as on social media:



.bock<sup>///</sup>  
**ACADEMY**





**.bock**<sup>®</sup> 

© registered trademark  
Updated: 07 | 2022  
Subject to technical changes.  
Picture credits: creanimation, Münster, Germany  
Nicola Rehage, Gütersloh

Hermann Bock GmbH  
Nickelstrasse 12  
33415 Verl  
Germany

Phone: +49 (0) 5246 9205 0  
Fax: +49 (0) 5246 9205 25  
Internet: [www.bock.net](http://www.bock.net)  
Email: [info@bock.net](mailto:info@bock.net)

domiflex<sup>®</sup> 3